

## Direct Flow Cytometry Protocol

Statement: The flow cytometry protocol provided is for reference only. A more suitable protocol should be designed based on the actual experimental context.

### Materials Required

- a. Samples: Cell lines or single-cell suspensions from tissues.
- b. Antibodies: Fluorochrome-conjugated primary antibodies (e.g., FITC, PE, APC), and isotype control antibodies.
- c. Buffers:
  - Staining buffer (Pre-cooled PBS + 1% BSA or 2% FBS).
  - Fixation buffer (1-4% paraformaldehyde or commercially available fixation reagents).
  - Permeabilization buffer (e.g., 0.1% Triton X-100 or commercial kits for intracellular staining).
- d. Other reagents: Fc receptor blocking solution (e.g., anti-CD16/32), viability dye (e.g., DAPI, 7-AAD).
- e. Equipment: Flow cytometer, centrifuge, refrigerator.
- f. Other materials: Centrifuge tubes/flow cytometry tubes, ice, pipettes.

### Procedure

#### 1. Planning Your Experiment

- a. Determine the specific cell populations and markers of interest.
- b. Select appropriate fluorochrome-conjugated antibodies that match with your flow cytometer's laser and filter setup.

#### 2. Cell Preparation

- a. Collect cells from culture or tissue as needed.
  - If isolating the cells from tissues, use an appropriate method such as enzymatic digestion and mechanical disruption.
  - Suspension cells: Collect the cell suspension into the centrifuge tube, centrifuge at  $300 \times g$  for 5 minutes, discard the supernatant and resuspend in staining buffer (remove the culture medium).
  - Adherent cells: Gently detach the cells from the culture flask or dish (trypsin-free cell dissociation solution is recommended), collect the cell suspension into the centrifuge tube, centrifuge at  $300 \times g$  for 5 minutes, discard the supernatant, and resuspend in staining buffer (remove the culture medium).

*Note: Please process the sample into a single-cell suspension.*

- b. Wash the cells with staining buffer, centrifuge at  $300 \times g$  for 5 minutes, and discard the supernatant. Repeat once.
- c. Count cell density and viability (with viability guaranteed to be at least 90%).
- d. Adjust cell density to  $1 \times 10^7$  cells/mL in staining buffer.
- e. Transfer 100  $\mu$ L of cell suspension (approximately  $1 \times 10^6$  cells) into flow cytometry tubes.

**3. Fc Receptor Blocking (Optional)**

If the samples you are testing are cells with high FcR expression, please proceed with this step.

- a. Add Fc blocking reagent to the cell suspension.
- b. Incubate 10-15 minutes at 4°C.

**4. Cell Staining for Surface Markers**

- a. Add the appropriate amount of antibody to each tube according to the manufacturer's instructions. Mix gently.
- b. Incubate for 30-60 minutes at 4°C in the dark.
- c. Wash the cells with staining buffer, centrifuge at  $300 \times g$  for 5 minutes at 4°C, discard the supernatant. Repeat once.
- d. Resuspend the cells in 200-500  $\mu$ L staining buffer for analysis.

*Note: If you are using an unlabeled primary antibody and a labeled secondary antibody, please refer to the indirect staining protocol.*

**5. Fixation (Optional)**

- a. Fix the cells with 100-500  $\mu$ L fixation buffer for 15 minutes at room temperature in the dark.
- b. Centrifuge and resuspend the cells in staining buffer for analysis.

**6. Intracellular Staining (If Required)**

- a. Fix the cells with 100-500  $\mu$ L fixation buffer for 15 minutes at room temperature in the dark.
- b. Wash the cells with staining buffer, centrifuge at  $300 \times g$  for 5 minutes, discard the supernatant. Repeat once.
- c. Permeabilization: Add 1-2 mL permeabilization buffer, mix gently and incubate at room temperature for 10-15 minutes in the dark.
- d. Wash the cells with permeabilization buffer, centrifuge at  $300 \times g$  for 5 minutes, discard the supernatant. Repeat once.
- e. Resuspend the cells in 100  $\mu$ L permeabilization buffer, add the appropriate amount of antibody to each tube according to the manufacturer's instructions. Mix gently and incubate for 30-60 minutes at 4°C in the dark.
- f. Wash the cells with permeabilization buffer, centrifuge at  $300 \times g$  for 5 minutes, discard the supernatant. Repeat once.
- g. Resuspend the cells in 200-500  $\mu$ L staining buffer for analysis.

**7. Viability Staining (Optional)**

Add viability dye (e.g., DAPI or 7-AAD) to the cells. Incubate for 10-15 minutes in the dark.

**8. Data Acquisition and Analysis**

- a. Prepare the flow cytometer by setting up voltage setting.
- b. Use unstained and single-stained controls for compensation.
- c. Collect  $\geq 10,000$  events per sample.
- d. Analyze with the flow cytometry software (e.g., FlowJo, FCS Express).