

Indirect Flow Cytometry Protocol

Statement: The flow cytometry protocol provided is for reference only. A more suitable protocol should be designed based on the actual experimental context.

Materials Required

- a. Samples: Cell lines or single-cell suspensions from tissues.
- b. Antibodies: Un-conjugated primary antibodies, fluorochrome-conjugated secondary antibodies (e.g., FITC, PE, APC), and isotype control antibodies.
- c. Buffers:
 - Staining buffer (Pre-cooled PBS + 1% BSA or 2% FBS).
 - Fixation buffer (1-4% paraformaldehyde or commercially available fixation reagents).
 - Permeabilization buffer (e.g., 0.1% Triton X-100 or commercial kits for intracellular staining).
- d. Other reagents: Fc receptor blocking solution (e.g., anti-CD16/32), viability dye (e.g., DAPI, 7-AAD).
- e. Equipment: Flow cytometer, centrifuge, refrigerator.
- f. Other materials: Centrifuge tubes/flow cytometry tubes, ice, pipettes.

Procedure

1. Planning Your Experiment

- a. Determine the specific cell populations and markers of interest.
- b. Select appropriate fluorochrome-conjugated antibodies that match with your flow cytometer's laser and filter setup.

2. Cell Preparation

- a. Collect cells from culture or tissue as needed.
 - If isolating the cells from tissues, use an appropriate method such as enzymatic digestion and mechanical disruption.
 - Suspension cells: Collect the cell suspension into the centrifuge tube, centrifuge at $300 \times g$ for 5 minutes, discard the supernatant and resuspend in staining buffer (remove the culture medium).
 - Adherent cells: Gently detach the cells from the culture flask or dish (trypsin-free cell dissociation solution is recommended), collect the cell suspension into the centrifuge tube, centrifuge at $300 \times g$ for 5 minutes, discard the supernatant, and resuspend in staining buffer (remove the culture medium).

Note: Please process the sample into a single-cell suspension.

- b. Wash the cells with staining buffer, centrifuge at $300 \times g$ for 5 minutes, and discard the supernatant.
- c. Count cell density and viability (with viability guaranteed to be at least 90%).
- d. Adjust cell density to 1×10^7 cells/mL in staining buffer.

Recombinant Proteins & Antibodies

- e. Transfer 100 μ L of cell suspension (approximately 1×10^6 cells) into flow cytometry tubes.

3. Fc Receptor Blocking (Optional)

If the samples you are testing are cells with high FcR expression, please proceed with this step.

- a. Add Fc blocking reagent to the cell suspension.
- b. Incubate 10-15 minutes at 4°C.

4. Cell Surface Staining

- a. Add the appropriate amount of the primary antibody to each tube according to the manufacturer's instructions. Mix gently.
- b. Incubate for 30-60 minutes at 4°C in the dark.
- c. Wash the cells with staining buffer, centrifuge at $300 \times g$ for 5 minutes at 4°C, discard the supernatant. Repeat once.
- d. Resuspend the cells in 100 μ L staining buffer, add the appropriate amount of the secondary antibody to each tube according to the manufacturer's instructions, mix gently. Incubate for 30-60 minutes at 4°C in the dark.
- e. Wash the cells with staining buffer, centrifuge at $300 \times g$ for 5 minutes at 4°C, discard the supernatant. Repeat once.
- f. Resuspend the cells in 200-500 μ L staining buffer for analysis.

5. Fixation (Optional)

- a. Fix the cells with 100-500 μ L fixation buffer for 15 minutes at room temperature in the dark.
- b. Centrifuge and resuspend the cells in staining buffer for analysis.

6. Intracellular Staining (If Required)

- a. Fix the cells with 100-500 μ L fixation buffer for 15 minutes at room temperature in the dark.
- b. Wash the cells with staining buffer, centrifuge at $300 \times g$ for 5 minutes, discard the supernatant. Repeat once.
- c. Permeabilization: Add 1-2 mL permeabilization buffer, mix gently and incubate at room temperature for 10-15 minutes in the dark.
- d. Wash the cells with permeabilization buffer, centrifuge at $300 \times g$ for 5 minutes, discard the supernatant. Repeat once.
- e. Resuspend the cells in 100 μ L permeabilization buffer, add the appropriate amount of the primary antibody to each tube according to the manufacturer's instructions. Mix gently and incubate for 30-60 minutes at 4°C in the dark.
- f. Wash the cells with permeabilization buffer, centrifuge at $300 \times g$ for 5 minutes, discard the supernatant. Repeat once.
- g. Resuspend the cells in 100 μ L permeabilization buffer, add the appropriate amount of the secondary antibody to each tube according to the manufacturer's instructions, mix gently. Incubate for 30-60 minutes at 4°C in the dark.



Recombinant Proteins & Antibodies

- h. Wash the cells with staining buffer, centrifuge at $300 \times g$ for 5 minutes at 4°C , discard the supernatant. Repeat once.
- i. Resuspend the cells in 200-500 μL staining buffer for analysis.

7. Viability Staining (Optional)

Add viability dye (e.g., DAPI or 7-AAD) to the cells. Incubate for 10-15 minutes in the dark.

8. Data Acquisition and Analysis

- a. Prepare the flow cytometer by setting up voltage setting.
- b. Use unstained and single-stained controls for compensation.
- c. Collect $\geq 10,000$ events per sample.
- d. Analyze with the flow cytometry software (e.g., FlowJo, FCS Express).